

免疫法大便檢體採集步驟

Stool blood testing of Immunology collect Procedure

1. 轉開綠色蓋子取出採便棒。

Screw open the green cap and take out the stirrer.

2. 在大便上不同的地方劃 6 下，少量的大便覆蓋溝槽即可。

Stroke 6 times on different spots of stool. Only a small amount of stool

3. 將綠色蓋子插回收集管，旋緊蓋子激烈搖盪混合。

Plug the green cap back into the collect tube. Screw it tight and shake fiercely to mix thoroughly.

■ 注意事項：Announcements

1. 在標籤寫上姓名、年齡及採便日期

Please write down your name, age, and collecting date on the label.

2. 採取大便檢體，採完後放入塑膠袋中，立即交付檢驗，否則需放在避光陰涼處。

Please look at the upper pictures step by step and collect stool. After collecting, put it in the plastic bag and take it to us in no time. But if you cannot do it right now, please keep it out of the light.

3. 請勿打開管瓶前端之錐形蓋。

Please don't open the front taper cover of the tube.

4. 請勿將管瓶中液體倒出，或任意在瓶中加水。

Please don't pour out the fluid in the tube or add water in it at will.

5. 大便硬時可加水弄濕，待大便軟化後再以採便棒插採，若有水便狀況，不易取樣時，請改日再取。

If the stool too hard to be collect, you can add some water to let it be soften. If you have aqueous stool and collect it difficultly, please collect it another day.

6. 遇痔瘡出血及女性月經期間，暫停採便。

If you have hemorrhoids with bleeding or in menstrual period, please don't collect stool.

7. 使用坐式馬桶時為避免大便掉入水中不易採取，請往前坐或反坐。並在便器內斜面先鋪上衛生紙以利採集。

When you use sitting stool, please sit forward or sit at the opposite way to avoid stool dropping into water. And put toilet paper on the ramp of stool for collecting well.

* 若有任何疑問，請電話詢問檢驗科

For any questions, please call laboratory for inquiry

正榮院區：(02) 24633330 轉分機 88445

Cheng-jueng District: (02)24633330#88445

孝二院區：(02) 24231735 轉分機 79117

Hsiao-er District: (02)24231735#79117

